

Flow Sensor Controlled Dispensing and Metering

Lübeck Summer Academy on Medical Technology

July 4, 2017 – Lübeck


Dr. Dorothee M. Runge, Dr. Thomas Weisener, Elisabeth Westphal

HNP Mikrosysteme GmbH • Bleicherufer 25 • 19053 Schwerin • T +49 385 52190-300 • info@hnp-mikrosysteme.de • www.hnp-mikrosysteme.de

Overview



- Controlled dispensing and metering
- Actuators
- Sensors
- Influencing factors on accuracy
- HNPM μ Dispense



- Actuator** Drive micro annular gear pump → revolution $n_{control}$
- Measuring element** Volume flow sensor
- Set value** Micro controller, set point setting by LabVIEW application (signal indicator)
- PID-controller** Implemented in micro controller

Accuracy, Trueness and Precision

According to DIN ISO 5725-2:2002-12 the term "accuracy" includes the two criteria "trueness" and "precision".

Trueness:

Measure for systematic deviation
How close are the results to the reference value?


$$R [\%] = \frac{\bar{V} - V_{Ref}}{V_{Ref}} \cdot 100$$

- R ... Trueness
- s ... Standard deviation
- V_{Ref} ... Reference value of volume
- \bar{V} ... Mean value of volume
- CV ... Coefficient of variation


Precision / Coefficient of Variation CV:

Scattering of the results
How close together are the results?


$$CV [\%] = \frac{s \cdot 100}{\bar{V}}$$




poor trueness
poor precision
very poor accuracy



good trueness
poor precision
poor accuracy



poor trueness
good precision
poor accuracy



good trueness
good precision
good accuracy

Dosing Pumps



Syringe pumps



Cavro



Cetoni



Hamilton

Piston pumps



FMI



Saphirwerk

HPLC-pumps



Knauer



PerkinElmer



Agilent

Annular gear pump



HNP Mikrosysteme

| pl | nl | μl | ml | l | % | ≈ | ⌚ | Δp | η | 📦 |
|----|----|----|----|---|---|---|---|----|---|---|
| | | | | | + | + | ○ | - | ○ | - |
| | | | | | + | + | ○ | ○ | ○ | - |
| | | | | | + | + | + | + | ○ | ○ |
| | | | | | + | + | + | ○ | + | + |

% precision ≈ resistance ⌚ life Δp pressure η Viscosity 📦 measurement
 low pulsation oscillating dispensing continuous delivery

Image source: named companies


Dosing Pumps



Drop-on-demand



Microdrop



GeSim

Si-diaphragm pumps



thinXS



Bartels



PARItec

Eccentric screw pump



ViscoTec

Osmotic pumps



Osmotex



Takasago

| pl | nl | μl | ml | l | % | ≈ | ⌚ | Δp | η | 📦 |
|----|----|----|----|---|---|---|---|----|---|---|
| | | | | | + | + | ○ | - | - | + |
| | | | | | ○ | ○ | - | - | - | + |
| | | | | | ○ | ○ | ○ | ○ | + | ○ |
| | | | | | - | - | - | ○ | + | + |

% precision ≈ resistance ⌚ life Δp pressure η viscosity 📦 measurement
 low pulsation oscillating dispensing continuous delivery

Image source: named companies

Volume Flow Sensors

Measuring Principles


Coriolis



Thermic



Differential Pressure



Ultrasound



nl / min μ l / min ml / min l / min

\square % \dot{m} f

| | | | | |
|--|--|--|--|--|
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |


\square measurement % precision \dot{m} mass flow (flow independent of liquid) f sampling rate

Image source: named companies

Dispense®

Applications in Life Science


Liquid Handling



Ion Channel Screening




Blood Typing



Sample Preparation




Flow Cytometry / Cell Sorting



In Vitro Diagnostics – IVD



DNA Isolation





Micro Dispense Module

- high precision of dispense volume and volume flow
- size of half-height syringe pump
- programming commands of syringe pump
- low noise level
- modular design

Micro annular gear pump

- long service life
- forward and reverse flow modes
- continuous volume flow
- wide volume flow range of one single pump

Sensor control and monitoring

- accuracy based on controlled volume flow or dosing volume
- fault detection (e.g. clogging, air bubbles)
- monitoring and documentation



Dispense

Dispense®

Modular Design




Sensor module
outlet (¼"-28UNF)
volume flow sensor
(up to 1000 µl/min)

Pump module
pump mZR-2521
(10 µl/min to 9 ml/min)

Valve module
2/2- way valve


3 Inlet modules
with filters
(¼"-28 UNF)




Micro Dispense Module

Cartridge Filter in a smallest format

- small size (L = 15 mm)
- low dead volume
- filter area 0,95 cm²
- mesh size 10 µm
- low pressure drop of 50 mbar at 200 ml/min with DI-water





| | |
|---|-------------|
| 1 | filter case |
| 2 | O-ring |
| 3 | filter |
| 4 | lid |
| 5 | washer |
| 6 | screw |
| 7 | fitting |
| 8 | tubing |




Dispense® Closed loop controlled volume flow

- **Closed loop controlled volume flow** by use of a flow sensor **10 µl/min to 1000 µl/min**
- Set value (T90) reached in max. 70 ms
- Precision 5% of measured value
- **Flow control without overshooting**



- Average dispensing precision μ Dispense with water
- Micro annular gear pump mZR-2521
 - Dispensing volume starting with 0,5 μ l
 - Volume flow 0,01 ml/min to 9 ml/min



Definition

The coefficient of variation (CV) expressed in % is defined as a standard deviation from the average value.

$$CV [\%] = \frac{s \cdot 100}{V}$$

CV Coefficient of variation [%]
 s Standard deviation [μ l]
 V Average volume [μ l]

Thank you for your attention

MECHANICAL
engineering

LIFE SCIENCE

CHEMICAL
engineering

» precise pumps – smart solutions «

micro annular gear pumps

HNP M
HNP Mikrosysteme

www.hnp-mikrosysteme.com

micro annular gear pumps

Pumping · Dosing · Metering

| | |
|-----------------------|---|
| Flow Rate | Q = 1 μ l/h ... 1 μ l/min ... 1 l/min |
| Dispense Volume | V = 0.1 ... 0.25 ... 1 ... 1000 μ l |
| Differential Pressure | Δp = -0.9 ... 40 ... 80 ... 150 bar |
| Precision | CV = 0.1 ... 0.5 ... 3 % |

mZR-2521